

The image features three medical catheters with inflatable balloons, arranged diagonally from the bottom left towards the top right. The catheters are made of clear plastic, and the balloons are inflated. The background is a light, neutral color with faint, curved lines that suggest a medical or scientific setting. The overall aesthetic is clean and professional.

Boa

Large diameter dilatation catheter

**Bigger Volume,
Wider Applications.**



Boa - Large diameter dilatation catheter.

Features.

- Available as up to a 20mm diameter balloon on a 6-7F low profile catheter.
- High pressure balloon inflation.
- Hydrophilic coating as standard.

Specifications.

Catheter Lengths: 75cm, 120cm, 150cm

Nominal pressure: 6atm

Rated Burst Pressure:	Balloon Diameter	Balloon Length: 20-50mm	Balloon Length: 60-80mm
	14	13atm	11atm
	16	13atm	11atm
	18	10atm	9atm
	20	9atm	8atm

Guide Wire Compatible: 0.035"

Coating: Hydrophilic

Product Codes.

Catheter Length (cm)	Catheter Diameter (Fr)	Introducer Size (Fr)	Balloon Diameter (mm)	Balloon Length (mm)					
				20	30	40	50	60	80
75	6	7	14	BOA142-75	BOA143-75	BOA144-75	BOA145-75	BOA146-75	BOA148-75
			16	BOA162-75	BOA163-75	BOA164-75	BOA165-75	BOA166-75	BOA168-75
	7	8	18	BOA182-75	BOA183-75	BOA184-75	BOA185-75	BOA186-75	BOA188-75
			20	BOA202-75	BOA203-75	BOA204-75	BOA205-75	BOA206-75	BOA208-75
120	6	7	14	BOA142-120	BOA143-120	BOA144-120	BOA145-120	BOA146-120	BOA148-120
			16	BOA162-120	BOA163-120	BOA164-120	BOA165-120	BOA166-120	BOA168-120
	7	8	18	BOA182-120	BOA183-120	BOA184-120	BOA185-120	BOA186-120	BOA188-120
			20	BOA202-120	BOA203-120	BOA204-120	BOA205-120	BOA206-120	BOA208-120
150	6	7	14	BOA142-150	BOA143-150	BOA144-150	BOA145-150	BOA146-150	BOA148-150
			16	BOA162-150	BOA163-150	BOA164-150	BOA165-150	BOA166-150	BOA168-150
	7	8	18	BOA182-150	BOA183-150	BOA184-150	BOA185-150	BOA186-150	BOA188-150
			20	BOA202-150	BOA203-150	BOA204-150	BOA205-150	BOA206-150	BOA208-150



Pan Medical Ltd

Barnett Way
Barnwood
Gloucester, UK
GL4 3RT
T: +44(0)1452 621 621
F: +44(0)1452 372 140
panmed@panmed.co.uk
www.panmed.co.uk