

# Clearpass

Self expandable peripheral stent



**Stent Longer.**



# Clearpass - Self expandable peripheral stent.

## Features.

- Minimal shortening after deployment.
- Radiopaque platinum iridium markers on the inner catheter.
- Controlled stent positioning.
- Extremely kink resistant and flexible.
- 3 integrated laser welded tantalum markers at each end.
- Two different working lengths 80cm / 120cm.
- Excellent pushability for cross-over.

## Specifications.

Guide Wire Compatible:

0.018", 0.021"

Recommended Sheath Compatibility:

5F

## Product Codes.

Catheter Length	Stent Length (mm)	Stent Diameter (mm)					
		5	6	7	8	9	10
80cm	30	CL53-80	CL63-80	CL73-80	-	-	-
	40	CL54-80	CL64-80	CL74-80	CL84-80	CL94-80	CL104-80
	50	CL55-80	CL65-80	CL75-80	CL85-80	CL95-80	CL105-80
	60	CL56-80	CL66-80	CL76-80	CL86-80	CL96-80	CL106-80
	80	CL58-80	CL68-80	CL78-80	CL88-80	CL98-80	CL108-80
	120	CL512-80	CL612-80	CL712-80	CL812-80	CL912-80	CL1012-80
	150	CL515-80	CL615-80	CL715-80	CL815-80	-	-
	180	CL518-80	CL618-80	CL718-80	-	-	-
	200	CL520-80	CL620-80	CL720-80	-	-	-
120cm	30	CL53-120	CL63-120	CL73-120	-	-	-
	40	CL54-120	CL64-120	CL74-120	CL84-120	CL94-120	CL104-120
	50	CL55-120	CL65-120	CL75-120	CL85-120	CL95-120	CL105-120
	60	CL56-120	CL66-120	CL76-120	CL86-120	CL96-120	CL106-120
	80	CL58-120	CL68-120	CL78-120	CL88-120	CL98-120	CL108-120
	120	CL512-120	CL612-120	CL712-120	CL812-120	CL912-120	CL1012-120
	150	CL515-120	CL615-120	CL715-120	CL815-120	-	-
	180	CL518-120	CL618-120	CL718-120	-	-	-
	200	CL520-120	CL620-120	CL720-120	-	-	-



### Pan Medical Ltd

Barnett Way  
Barnwood  
Gloucester, UK  
GL4 3RT

T: +44(0)1452 621 621  
F: +44(0)1452 372 140  
panmed@panmed.co.uk  
www.panmed.co.uk