

Courageous

Extra long dilatation catheter

**Go Deeper.
Go Longer.**



Courageous - Extra Long RX and OTW Dilatation Catheters.

Features.

- For all 0.014" and 0.018" peripheral needs.
- Available as both Over The Wire (OTW) and Rapid Exchange (RX) catheters.
- Up to 300mm balloon length and 22atm high pressure balloons.

Specifications.

Average Compliance:	7%												
Catheter Lengths:	100cm, 130cm, 150cm for OTW — 140cm for RX												
Nominal pressure:	6atm												
Rated Burst Pressure:													
	<table border="1"> <thead> <tr> <th>Balloon Diameter</th> <th>Balloon Length: 100-150mm</th> <th>Balloon Length: 200-250mm</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>22atm</td> <td>20atm</td> </tr> <tr> <td>4</td> <td>22atm</td> <td>20atm</td> </tr> <tr> <td>5</td> <td>22atm</td> <td>20atm</td> </tr> </tbody> </table>	Balloon Diameter	Balloon Length: 100-150mm	Balloon Length: 200-250mm	3	22atm	20atm	4	22atm	20atm	5	22atm	20atm
Balloon Diameter	Balloon Length: 100-150mm	Balloon Length: 200-250mm											
3	22atm	20atm											
4	22atm	20atm											
5	22atm	20atm											
Guide Wire Compatible:	0.014" and 0.018"												
Sheath Compatibility:	4F (OTW only)												

Product Codes.

Catheter Length	Balloon Diameter (mm)	Balloon Length (mm)						
		100	120	150	200	250	300	
100cm	3	CR310-100	CR312-100	CR315-100	CR320-100	CR325-100	CR330-100	OTW
100cm	4	CR410-100	CR412-100	CR415-100	CR420-100	CR425-100	CR430-100	
100cm	5	CR510-100	CR512-100	CR515-100	CR520-100	CR525-100	CR530-100	
130cm	3	CR310-130	CR312-130	CR315-130	CR320-130	CR325-130	CR330-130	
130cm	4	CR410-130	CR412-130	CR415-130	CR420-130	CR425-130	CR430-130	
130cm	5	CR510-130	CR512-130	CR515-130	CR520-130	CR525-130	CR530-130	
150cm	3	CR310-150	CR312-150	CR315-150	CR320-150	CR325-150	CR330-150	
150cm	4	CR410-150	CR412-150	CR415-150	CR420-150	CR425-150	CR430-150	
150cm	5	CR510-150	CR512-150	CR515-150	CR520-150	CR525-150	CR530-150	
140cm	3	CRRX310-140	CRRX312-140	CRRX315-140	CRRX320-140	CRRX325-140	CRRX330-140	RX
140cm	4	CRRX410-140	CRRX412-140	CRRX415-140	CRRX420-140	CRRX425-140	CRRX430-140	
140cm	5	CRRX510-140	CRRX512-140	CRRX515-140	CRRX520-140	CRRX525-140	CRRX530-140	

