



PERIBALT

Cobalt Chromium Balloon Expandable Peripheral Stent



PERIBALT

Cobalt Chromium Balloon Expandable Peripheral Stent

The Peribalt is a pre mounted peripheral stent made from Cobalt Chromium.

It offers all the advantages of Cobalt Chromium in Peripheral Stenting.

Features.

- The Peribalt has been developed with up to a 35% reduction in strut thickness.
- The cobalt chromium alloys are far more radiopaque than stainless steel - for precise visualisation and positioning.
- The stent structure has a low profile - providing improved flexibility.
- Pre-mounted on a tri-folded balloon – to allow uniform stent expansion.
- Low profile of crimped stent – allows outstanding delivery capabilities.

Excellent radial strength.

Minimum shortening of expanded stent.

Ordering Information.

Stent Length (mm)	Stent Diameter (mm)					
	5	6	7	8	9	10
12	PA/PB/12/05	PA/PB/12/06	PA/PB/12/07	PA/PB/12/08	PA/PB/12/09	PA/PB/12/10
16	PA/PB/16/05	PA/PB/16/06	PA/PB/16/07	PA/PB/16/08	PA/PB/16/09	PA/PB/16/10
25	PA/PB/25/05	PA/PB/25/06	PA/PB/25/07	PA/PB/25/08	PA/PB/25/09	PA/PB/25/10
35	PA/PB/35/05	PA/PB/35/06	PA/PB/35/07	PA/PB/35/08	PA/PB/35/09	PA/PB/35/10
45	PA/PB/45/05	PA/PB/45/06	PA/PB/45/07	PA/PB/45/08	PA/PB/45/09	PA/PB/45/10
55	PA/PB/55/05	PA/PB/55/06	PA/PB/55/07	PA/PB/55/08	PA/PB/55/09	PA/PB/55/10
72	PA/PB/72/05	PA/PB/72/06	PA/PB/72/07	PA/PB/72/08	PA/PB/72/09	PA/PB/72/10

All stents are mounted on a 75cm catheter with a 5F shaft



PAN Medical
Interventional Radiology

Barnett Way
Barnwood
Gloucester
GL4 3RT, UK

Tel: +44 (0)1452 621621
Fax: +44 (0)1452 372140
Email: panmed@panmed.co.uk
Website: www.panmed.co.uk