

PHANTOM GLIDE

Fully coated PTA Catheter



Features:

- Designed for challenging PTA procedures
- High RBP
- Excellent re-wrapping

Ordering Information:

Specifications:

Balloon Diameters:	3 - 14mm
Balloon Lengths:	20 - 120mm
Catheter Lengths:	50, 75, 135cm
Sheath Compatibility:	5 - 7F
Guide Wire Compatibility:	"0.035"
RBP:	Up to 22ATM
Balloon Material:	Nylon Polymer (average compliance 7%)
Coating:	Hydrophilic
Balloon Fold:	Quint Fold (5 folds)

Shaft Length:			Nominal Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (Fr)
50cm	75cm	135cm			
PG-3250	PG-3275	PG-32135	3	20	5
PG-4250	PG-4275	PG-42135	4	20	5
PG-5250	PG-5275	PG-52135	5	20	5
PG-6250	PG-6275	PG-62135	6	20	5
PG-7250	PG-7275	PG-72135	7	20	5
PG-8250	PG-8275	PG-82135	8	20	6
PG-3350	PG-3375	PG-33135	3	30	5
PG-4350	PG-4375	PG-43135	4	30	5
PG-5350	PG-5375	PG-53135	5	30	5
PG-6350	PG-6375	PG-63135	6	30	5
PG-7350	PG-7375	PG-73135	7	30	5
PG-8350	PG-8375	PG-83135	8	30	6
PG-9350	PG-9375	PG-93135	9	30	6
PG-10350	PG-10375	PG-103135	10	30	7
PG-3450	PG-3475	PG-34135	3	40	5
PG-4450	PG-4475	PG-44135	4	40	5
PG-5450	PG-5475	PG-54135	5	40	5
PG-6450	PG-6475	PG-64135	6	40	5
PG-7450	PG-7475	PG-74135	7	40	5
PG-8450	PG-8475	PG-84135	8	40	6
PG-9450	PG-9475	PG-94135	9	40	6
PG-10450	PG-10475	PG-104135	10	40	7
PG-12450	PG-12475	PG-124135	12	40	7
PG-3550	PG-3575	PG-35135	3	50	5
PG-4550	PG-4575	PG-45135	4	50	5
PG-5550	PG-5575	PG-55135	5	50	5
PG-6550	PG-6575	PG-65135	6	50	5
PG-7550	PG-7575	PG-75135	7	50	5
PG-8550	PG-8575	PG-85135	8	50	6
PG-9550	PG-9575	PG-95135	9	50	6
PG-10550	PG-10575	PG-105135	10	50	7
PG-12550	PG-12575	PG-125135	12	50	7
PG-3650	PG-3675	PG-36135	3	60	5
PG-4650	PG-4675	PG-46135	4	60	5
PG-5650	PG-5675	PG-56135	5	60	5
PG-6650	PG-6675	PG-66135	6	60	5
PG-7650	PG-7675	PG-76135	7	60	5
PG-8650	PG-8675	PG-86135	8	60	6
PG-9650	PG-9675	PG-96135	9	60	6
PG-10650	PG-10675	PG-106135	10	60	7
PG-12650	PG-12675	PG-126135	12	60	7
PG-14650	PG-14675	PG-146135	14	60	7
PG-3850	PG-3875	PG-38135	3	80	5
PG-4850	PG-4875	PG-48135	4	80	5
PG-5850	PG-5875	PG-58135	5	80	5
PG-6850	PG-6875	PG-68135	6	80	5
PG-7850	PG-7875	PG-78135	7	80	5
PG-8850	PG-8875	PG-88135	8	80	6
PG-9850	PG-9875	PG-98135	9	80	6
PG-10850	PG-10875	PG-108135	10	80	7
PG-12850	PG-12875	PG-128135	12	80	7
PG-14850	PG-14875	PG-148135	14	80	7
PG-31050	PG-31075	PG-310135	3	100	5
PG-41050	PG-41075	PG-410135	4	100	5
PG-51050	PG-51075	PG-510135	5	100	5
PG-61050	PG-61075	PG-610135	6	100	5
PG-71050	PG-71075	PG-710135	7	100	5
PG-81050	PG-81075	PG-810135	8	100	6
PG-91050	PG-91075	PG-910135	9	100	6
PG-101050	PG-101075	PG-1010135	10	100	7
PG-121050	PG-121075	PG-1210135	12	100	7
PG-141050	PG-141075	PG-1410135	14	100	7
PG-31250	PG-31275	PG-312135	3	120	5
PG-41250	PG-41275	PG-412135	4	120	5
PG-51250	PG-51275	PG-512135	5	120	5
PG-61250	PG-61275	PG-612135	6	120	5
PG-71250	PG-71275	PG-712135	7	120	5
PG-81250	PG-81275	PG-812135	8	120	6
PG-91250	PG-91275	PG-912135	9	120	6
PG-101250	PG-101275	PG-1012135	10	120	7
PG-121250	PG-121275	PG-1212135	12	120	7
PG-141250	PG-141275	PG-1412135	14	120	7

Other sizes are available upon request



PAN Medical
Interventional Radiology

Barnett Way
Barnwood
Gloucester
GL4 3RT, UK

Tel: +44 (0)1452 621621
Fax: +44 (0)1452 372140
Email: panmed@panmed.co.uk
Website: www.panmed.co.uk